

Two sets of capture probes were prepared for detection of cucumber mosaic virus.

CMV capture probe sequences

Probe set #1: Recognizes region nt 1542 – 1641 in the CMV genome	
3' probe: CMV-8c.1	CCCACACGGTAGAATCAAATTTTCGGCAAAGGATTA
5' probe: CMV-8c.2	ACTCGAATTTGAATGCGCGAAACAAGCTTCTTATC
Probe set #2: Recognizes region nt 140 – 271 in the CMV genome	
3' probe: CMV-11c.2	CAATGATGCCGGCCTAGGTCACACTCAGTAGCCAT
5' probe: CMV-6c.2	TTTGAAGATCGTCAGACGTAGCCGCTGAGGACTGT

Sequences from Agindotan and Perry, *Phytochemistry*, **97** 1 (2007).